

Tips for HPSA's Medical Sharps Collection Program

Tip #1

Make sure the Fill Line is on the same side as the BIOHAZARD diamond.

Since the medical sharps container box is rectangular it is easy to make the mistake of putting the lid on backwards. When installed correctly, **the Fill Line should be on the same side as the BIOHAZARD diamond.** This ensures that the container handle is not an obstruction to the user when inserting medical sharps.



Tip #2

Make sure to secure the lid to the base.

Remember to secure the lid to base of the container, by **pushing down on each corner of the lid until you hear a click.**



Tip #3

Make sure to temporarily close your container.

Lightly press the container lid shut until you hear a click. **The rectangular tab will remain in an upright position for you to easily reopen your container.**



Tip #4

Make sure to permanently close your container.

Press down firmly on the 2 textured divots located on the center of the lid until you hear a click. **When closed permanently, the rectangular tab should no longer be accessible to open the container.**



Reference:
[MedicalSharpsTips/Allprovinces/2024/V1](#)

Need more help? Refer to our [video](#) to view how to assemble the Medical Sharps Container.



For any questions, please contact us at:
info@healthsteward.ca
Visit us at: healthsteward.ca